

How to use the safety straps:


A-Frame - Standard Crutches: Straps are intended only as additional support when carrying heavier items. Bring straps around the hand grip from front of the crutch to the back. Loop and hook strips are located in such a way to allow safe looping around hand grip with minimal interference of a proper grip. Attach the ends of the straps to hook portion location inside the elastic backing. They are located in such a way that the frame of the crutch will press against them for an even stronger grip and added security.

Walkers, Wheelchairs and other items: Loop each strap up over and around the pole so that the strap loops back up onto itself and two Velcro sides meet at the back of the bag. Once straps have been attached to backing, tuck it into the elastic backing for added security.

How to use your new Krutch Kaddy with standard A Frame crutches:



1. Slip the foot of the crutch through the elastic backing located on the rear of your new Krutch Kaddy.



3. As you get closer to the hand grips, stretch out the elastic backing, this will make it easier to slide the Krutch Kaddy into position.



2. Grasp the top of both the front and the elastic backing. Slide the Krutch Kaddy up along the crutch frame.



4. Position the Krutch Kaddy where it is most comfortable for you. We recommend 1-2" below hand grip.

How to care for your Krutch Kaddy:

You can machine-wash separately in cold water using a mild detergent.

Cleaning: Your Krutch Kaddy can be machine washed in the gentle cycle. Do not bleach. To maintain the shape we recommend air drying only.

Not just for crutches: the Krutch Kaddy® can be used with any item that has a pole on which you can attach the two Velcro straps

How to use with forearm crutches:

- If you are able to easily remove the cuff attachment:



1. Turn the Plastic Locking Nut counter clockwise until it slide freely up the pole. Push down the buttons located on each side of the crutch pole.



4. Reattach the cuff to the crutch pole making sure the buttons lock into place before using forearm crutches.



2. Grasp the cuff and pull up off of crutch pole.



5. Extend the two straps located on the back of the Krutch Kaddy. Loop each strap up over and around the cuff so that the strap loops back up onto itself and two Velcro sides meet at the back of the bag.



3. Slide the portion of the crutch pole located above the hand grip up through the elastic backing on the back of your Krutch Kaddy.



6. You may have some "extra" strap showing, you can tuck it into the elastic backing as shown. Your Krutch Kaddy is now securely in place and ready

NOTE: If you are unable to easily remove the cuff attachment: Slide the foot of the crutch down between the elastic and the back of the Kaddy. Slide the bag up the pole. Getting the bag over the hand grip will take some effort. I needed two people. Stretch the elastic so that you can get the bag over the hand grip and onto the pole section between the cuff and the grip. Then loop straps up over the cuff and secure at back of Kaddy.