

# Push<sup>®</sup> med Wrist Brace

Immobilizes wrist with less bulk, greater comfort

- Anatomically shaped dorsal and volar stays fit the exact contours of your hand and wrist
- Wide elastic band permits precise adjustment
- Design avoids pressure to vulnerable bony wrist areas
- Brace is easily applied and removed with one hand
- Streamlined contour allows for easy dressing and undressing
- Can be machine washed on gentle cycle at low temperature
- Creates soft, comfortable feel and wicks moisture through the Sympress<sup>™</sup> microfibers



## Indications

- Sprains / strains or instability of the wrist
- Tendinitis of wrist tendon(s)
- Non-displaced fracture
- Support after wrist fracture
- Arthritis of the wrist

## Measuring / Sizing

Measure around the wrist as shown and choose a size from the chart. Note: If you have a cone shaped forearm you may need to choose a larger size.



### Sizing: Right or Left Wrist

SIZE	INCHES	CENTIMETERS
1	5 1/8 - 6	13 - 15
2	6 - 6 5/8	15 - 17
3	6 5/8 - 7 1/2	17 - 19
4	7 1/2 - 8 1/4	19 - 21



*The Push med Wrist Brace must be correctly SIZED and WRAPPED to work! Visit [BraceLab.com](http://BraceLab.com) for more information on fitting and care or consult your therapist.*

## Insurance Information

- The Push med Wrist Brace has an assigned L-code of L-3908.
- BraceLab is not a provider. We do not file insurance nor can we provide information on insurance coverage or reimbursement rates. Please contact your insurance for filing information.



Available in the US from [www.bracelab.com](http://www.bracelab.com)  
[www.bracelab.com](http://www.bracelab.com) · P: 601-892-3115 · F: 601-892-3116

