


**THE WRIGHT STUFF**




Health Care  
Products  
that make life a  
little easier

## ELECTRIC-RAZOR CUFF



1. Slide cuff over hand. Use Cuff Finger Loop to adjust strap to desired fit around palm of hand.
2. Insert electric razor into front Strap Pouch. Use front Finger Loops adjust straps to hold razor securely. Top Strap may also be used depending on the razor style.



08819

 **Therafin**  
Corporation

9450 W Laraway Rd • Frankfort, IL 60423  
Customer 800-843-7234 • Fax 888-479-1515  
[www.Therafin.com](http://www.Therafin.com) • [info@Therafin.com](mailto:info@Therafin.com)